

Array Layouts: Human Genome RT² miRNA PCR Array, 96-Well

MAH-001

001	1	2	3	4	5	6	7	8	9	10	11	12
A	miR-142-5p	miR-16	miR-142-3p	miR-21	miR-15a	miR-29b	let-7a	miR-126	miR-143	let-7b	miR-27a	let-7f
B	miR-9	miR-26a	miR-24	miR-30e	miR-181a	miR-29a	miR-124	miR-144	miR-30d	miR-19b	miR-22	miR-122
C	miR-150	miR-32	miR-155	miR-140-5p	miR-125b	miR-141	miR-92a	miR-424	miR-191	miR-17	miR-130a	miR-20a
D	miR-27b	miR-26b	miR-146a	miR-200c	miR-99a	miR-19a	miR-23a	miR-30a	let-7i	miR-93	let-7c	miR-106b
E	miR-101	let-7g	miR-425	miR-15b	miR-28-5p	miR-18a	miR-25	miR-23b	miR-302a	miR-186	miR-29c	miR-7
F	let-7d	miR-30c	miR-181b	miR-223	miR-320	miR-374a	let-7e	miR-151-5p	miR-374b	miR-196b	miR-140-3p	miR-100
G	miR-103	miR-96	miR-302b	miR-194	miR-125a-5p	miR-423-5p	miR-376c	miR-195	miR-222	miR-28-3p	miR-128a	miR-302c
H	miR-423-3p	miR-185	miR-30b	miR-210	SNORD48	SNORD47	SNORD44	U6	miRTC	miRTC	PPC	PPC

MAH-002

002	1	2	3	4	5	6	7	8	9	10	11	12
A	mir-375	mir-182	mir-196a	mir-10a	mir-324-5p	mir-137	mir-378	mir-135b	mir-342-3p	mir-205	mir-192	mir-10b
B	mir-744	mir-363	mir-503	mir-130b	mir-151-3p	mir-214	mir-454	mir-379	mir-301a	mir-98	mir-34a	mir-33a
C	mir-339-3p	mir-301b	mir-193b	mir-652	mir-20b	mir-138	mir-199a-3p	mir-183	mir-134	mir-132	mir-148a	mir-218
D	mir-365	mir-339-5p	mir-497	mir-193a-5p	mir-625	mir-488	mir-129-5p	mir-345	mir-518b	mir-517a	mir-421	mir-197
E	mir-17*	mir-18b	mir-335	mir-376a	mir-361-5p	mir-505	mir-877	mir-203	mir-411	mir-152	mir-532-5p	mir-342-5p
F	mir-373	mir-135a	mir-202	mir-340	mir-139-5p	mir-129-3p	mir-545	mir-338-3p	mir-127-3p	mir-455-3p	mir-484	mir-500
G	mir-518c	mir-1	mir-590-3p	mir-372	mir-148b	mir-181c	mir-188-5p	mir-149	mir-410	mir-431	mir-331-3p	mir-548c-5p
H	mir-486-5p	mir-34c-5p	mir-629	mir-452	SNORD48	SNORD47	SNORD44	U6	miRTC	miRTC	PPC	PPC

MAH-003

003	1	2	3	4	5	6	7	8	9	10	11	12
A	mir-574-5p	mir-184	mir-660	mir-337-5p	mir-146b-5p	mir-432	mir-370	mir-95	mir-206	mir-873	mir-532-3p	mir-139-3p
B	mir-524-5p	mir-188-3p	mir-590-5p	mir-455-5p	mir-331-5p	mir-429	mir-362-5p	mir-525-3p	mir-493	mir-526b	mir-190	mir-515-5p
C	mir-665	mir-147	mir-361-3p	mir-514	mir-376b	mir-582-5p	mir-187	mir-671-5p	mir-487b	mir-651	mir-219-5p	mir-299-3p
D	mir-449a	mir-369-3p	mir-219-2-3p	mir-105	mir-551b	mir-542-5p	mir-542-3p	mir-618	mir-502-3p	mir-519a	mir-181d	mir-382
E	mir-508-3p	mir-450a	mir-127-5p	mir-523	mir-589	mir-369-5p	mir-615-3p	mir-527	mir-576-3p	mir-337-3p	mir-362-3p	mir-515-3p
F	mir-518a-3p	mir-885-5p	mir-433	mir-512-3p	mir-548d-5p	mir-133b	mir-296-5p	mir-371-3p	mir-330-3p	mir-329	mir-495	mir-450b-5p
G	mir-758	mir-522	mir-323-3p	mir-889	mir-208b	mir-487a	mir-628-3p	mir-628-5p	mir-671-3p	mir-485-3p	mir-548d-3p	mir-383
H	mir-548b-3p	mir-215	mir-520d-5p	mir-520d-3p	SNORD48	SNORD47	SNORD44	U6	miRTC	miRTC	PPC	PPC

MAH-004

004	1	2	3	4	5	6	7	8	9	10	11	12
A	mir-584	mir-513-5p	mir-579	mir-498	mir-550	mir-556-5p	mir-380	mir-491-3p	mir-576-5p	mir-597	mir-504	mir-520g
B	mir-642	mir-875-5p	mir-875-3p	mir-501-5p	mir-501-3p	mir-544	mir-548a-5p	mir-548a-3p	mir-422a	mir-489	mir-636	mir-198
C	mir-496	mir-518d-3p	mir-524-3p	mir-519b-5p	mir-519c-3p	mir-516a-3p	mir-502-5p	mir-512-5p	mir-566	mir-548c-3p	mir-555	mir-492
D	mir-638	mir-412	mir-588	mir-650	mir-619	mir-543	mir-643	mir-603	mir-602	mir-563	mir-558	mir-657
E	mir-511	mir-595	mir-613	mir-802	mir-583	mir-384	mir-611	mir-575	mir-658	mir-637	mir-580	mir-564
F	mir-610	mir-633	mir-562	mir-614	mir-325	mir-600	mir-630	mir-298	mir-622	mir-626	mir-659	mir-585
G	mir-639	mir-220c	mir-211	mir-623	mir-641	mir-557	mir-212	mir-632	mir-567	mir-605	mir-208	mir-448
H	mir-621	mir-648	mir-190b	mir-587	SNORD48	SNORD47	SNORD44	U6	miRTC	miRTC	PPC	PPC

Gene Tables: MAH-001

Position	mirBASE v10.0 Accession Number	miRNA ID	RT ² miRNA qPCR Assay Cat. #	Position	mirBASE v10.0 Accession Number	miRNA ID	RT ² miRNA qPCR Assay Cat. #
A01	MIMAT0000433	hsa-miR-142-5p	MPH00043A	E01	MIMAT0000099	hsa-miR-101	MPH00011A
A02	MIMAT0000069	hsa-miR-16	MPH00062A	E02	MIMAT0000414	hsa-let-7g	MPH00007A
A03	MIMAT0000434	hsa-miR-142-3p	MPH01176A	E03	MIMAT0003393	hsa-miR-425	MPH00212A
A04	MIMAT0000076	hsa-miR-21	MPH00106A	E04	MIMAT0000417	hsa-miR-15b	MPH00061A
A05	MIMAT0000068	hsa-miR-15a	MPH00060A	E05	MIMAT0000085	hsa-miR-28-5p	MPH00134A
A06	MIMAT0000100	hsa-miR-29b	MPH01245A	E06	MIMAT0000072	hsa-miR-18a	MPH00075A
A07	MIMAT0000062	hsa-let-7a	MPH00001A	E07	MIMAT0000081	hsa-miR-25	MPH01235A
A08	MIMAT0000445	hsa-miR-126	MPH01161A	E08	MIMAT0000418	hsa-miR-23b	MPH01233A
A09	MIMAT0000435	hsa-miR-143	MPH01177A	E09	MIMAT0000684	hsa-miR-302a	MPH01247A
A10	MIMAT0000063	hsa-let-7b	MPH00002A	E10	MIMAT0000456	hsa-miR-186	MPH00072A
A11	MIMAT0000084	hsa-miR-27a	MPH01239A	E11	MIMAT0000681	hsa-miR-29c	MPH01246A
A12	MIMAT0000067	hsa-let-7f	MPH00006A	E12	MIMAT0000252	hsa-miR-7	MPH00427A
B01	MIMAT0000441	hsa-miR-9	MPH00456A	F01	MIMAT0000065	hsa-let-7d	MPH00004A
B02	MIMAT0000082	hsa-miR-26a	MPH00130A	F02	MIMAT0000244	hsa-miR-30c	MPH00152A
B03	MIMAT0000080	hsa-miR-24	MPH01234A	F03	MIMAT0000257	hsa-miR-181b	MPH00065A
B04	MIMAT0000692	hsa-miR-30e	MPH00154A	F04	MIMAT0000280	hsa-miR-223	MPH01231A
B05	MIMAT0000256	hsa-miR-181a	MPH00064A	F05	MIMAT0000510	hsa-miR-320	MPH00157A
B06	MIMAT0000086	hsa-miR-29a	MPH01244A	F06	MIMAT0000727	hsa-miR-374a	MPH00189A
B07	MIMAT0000422	hsa-miR-124	MPH01157A	F07	MIMAT0000066	hsa-let-7e	MPH00005A
B08	MIMAT0000436	hsa-miR-144	MPH01178A	F08	MIMAT0004697	hsa-miR-151-5p	MPH00055A
B09	MIMAT0000245	hsa-miR-30d	MPH00153A	F09	MIMAT0004955	hsa-miR-374b	MPH00190A
B10	MIMAT0000074	hsa-miR-19b	MPH01215A	F10	MIMAT0001080	hsa-miR-196b	MPH00086A
B11	MIMAT0000077	hsa-miR-22	MPH01228A	F11	MIMAT0004597	hsa-miR-140-3p	MPH01174A
B12	MIMAT0000421	hsa-miR-122	MPH00020A	F12	MIMAT0000098	hsa-miR-100	MPH00009A
C01	MIMAT0000451	hsa-miR-150	MPH00054A	G01	MIMAT0000101	hsa-miR-103	MPH00012A
C02	MIMAT0000090	hsa-miR-32	MPH00156A	G02	MIMAT0000095	hsa-miR-96	MPH00479A
C03	MIMAT0000646	hsa-miR-155	MPH00059A	G03	MIMAT0000715	hsa-miR-302b	MPH01248A
C04	MIMAT0000431	hsa-miR-140-5p	MPH00041A	G04	MIMAT0000460	hsa-miR-194	MPH00083A
C05	MIMAT0000423	hsa-miR-125b	MPH00023A	G05	MIMAT0000443	hsa-miR-125a-5p	MPH00022A
C06	MIMAT0000432	hsa-miR-141	MPH01175A	G06	MIMAT0004748	hsa-miR-423-5p	MPH00210A
C07	MIMAT0000092	hsa-miR-92a	MPH01375A	G07	MIMAT0000720	hsa-miR-376c	MPH00195A
C08	MIMAT0001341	hsa-miR-424	MPH00211A	G08	MIMAT0000461	hsa-miR-195	MPH00084A
C09	MIMAT0000440	hsa-miR-191	MPH00079A	G09	MIMAT0000279	hsa-miR-222	MPH01230A
C10	MIMAT0000070	hsa-miR-17	MPH00063A	G10	MIMAT0004502	hsa-miR-28-3p	MPH01241A
C11	MIMAT0000425	hsa-miR-130a	MPH01165A	G11	MIMAT0000424	hsa-miR-128a	MPH00026A
C12	MIMAT0000075	hsa-miR-20a	MPH00104A	G12	MIMAT0000717	hsa-miR-302c	MPH01249A
D01	MIMAT0000419	hsa-miR-27b	MPH01240A	H01	MIMAT0001340	hsa-miR-423-3p	MPH01289A
D02	MIMAT0000083	hsa-miR-26b	MPH00131A	H02	MIMAT0000455	hsa-miR-185	MPH00071A
D03	MIMAT0000449	hsa-miR-146a	MPH00047A	H03	MIMAT0000420	hsa-miR-30b	MPH00151A
D04	MIMAT0000617	hsa-miR-200c	MPH01218A	H04	MIMAT0000267	hsa-miR-210	MPH00107A
D05	MIMAT0000097	hsa-miR-99a	MPH00481A	H05	N/A	hsa-SNORD48	MPH01657A
D06	MIMAT0000073	hsa-miR-19a	MPH01214A	H06	N/A	hsa-SNORD47	MPH01660A
D07	MIMAT0000078	hsa-miR-23a	MPH01232A	H07	N/A	hsa-SNORD44	MPH01658A
D08	MIMAT0000087	hsa-miR-30a	MPH00150A	H08	N/A	hsa-U6	MPH01653A
D09	MIMAT0000415	hsa-let-7i	MPH00008A	H09	N/A	miRTC	N/A
D10	MIMAT0000093	hsa-miR-93	MPH00465A	H10	N/A	miRTC	N/A
D11	MIMAT0000064	hsa-let-7c	MPH00003A	H11	N/A	PPC	N/A
D12	MIMAT0000680	hsa-miR-106b	MPH00015A	H12	N/A	PPC	N/A

Gene Table: MAH-002

Position	mirBASE v10.0 Accession Number	miRNA ID	RT ² miRNA qPCR Assay Cat. #	Position	mirBASE v10.0 Accession Number	miRNA ID	RT ² miRNA qPCR Assay Cat. #
A01	MIMAT0000728	hsa-miR-375	MPH00191A	E01	MIMAT0000071	hsa-mir-17*	MPH01193A
A02	MIMAT0000259	hsa-miR-182	MPH00068A	E02	MIMAT0001412	hsa-miR-18b	MPH00076A
A03	MIMAT0000226	hsa-miR-196a	MPH00085A	E03	MIMAT0000765	hsa-miR-335	MPH00166A
A04	MIMAT0000253	hsa-miR-10a	MPH00017A	E04	MIMAT0000729	hsa-miR-376a	MPH00193A
A05	MIMAT0000761	hsa-miR-324-5p	MPH00159A	E05	MIMAT0000703	hsa-miR-361-5p	MPH00179A
A06	MIMAT0000429	hsa-miR-137	MPH00038A	E06	MIMAT0002876	hsa-miR-505	MPH01310A
A07	MIMAT0000732	hsa-miR-378	MPH01284A	E07	MIMAT0004949	hsa-miR-877	MPH00445A
A08	MIMAT0000758	hsa-miR-135b	MPH00036A	E08	MIMAT0000264	hsa-miR-203	MPH00098A
A09	MIMAT0000753	hsa-miR-342-3p	MPH01270A	E09	MIMAT0003329	hsa-miR-411	MPH00206A
A10	MIMAT0000266	hsa-miR-205	MPH00100A	E10	MIMAT0000438	hsa-miR-152	MPH00056A
A11	MIMAT0000222	hsa-miR-192	MPH00080A	E11	MIMAT0002888	hsa-miR-532-5p	MPH00293A
A12	MIMAT0000254	hsa-miR-10b	MPH00018A	E12	MIMAT0004694	hsa-miR-342-5p	MPH00173A
B01	MIMAT0004945	hsa-miR-744	MPH00428A	F01	MIMAT0000726	hsa-miR-373	MPH01280A
B02	MIMAT0000707	hsa-miR-363	MPH01276A	F02	MIMAT0000428	hsa-miR-135a	MPH00035A
B03	MIMAT0002874	hsa-miR-503	MPH00248A	F03	MIMAT0002811	hsa-miR-202	MPH01219A
B04	MIMAT0000691	hsa-miR-130b	MPH01166A	F04	MIMAT0004692	hsa-miR-340	MPH00172A
B05	MIMAT0000757	hsa-miR-151-3p	MPH01186A	F05	MIMAT0000250	hsa-miR-139-5p	MPH00040A
B06	MIMAT0000271	hsa-miR-214	MPH01223A	F06	MIMAT0004605	hsa-miR-129-3p	MPH01164A
B07	MIMAT0003885	hsa-miR-454	MPH01296A	F07	MIMAT0003165	hsa-miR-545	MPH01340A
B08	MIMAT0000733	hsa-miR-379	MPH00198A	F08	MIMAT0000763	hsa-miR-338-3p	MPH01265A
B09	MIMAT0000688	hsa-miR-301a	MPH00144A	F09	MIMAT0000446	hsa-miR-127-3p	MPH01162A
B10	MIMAT0000096	hsa-miR-98	MPH00480A	F10	MIMAT0004784	hsa-miR-455-3p	MPH01297A
B11	MIMAT0000255	hsa-miR-34a	MPH00176A	F11	MIMAT0002174	hsa-miR-484	MPH00228A
B12	MIMAT0000091	hsa-miR-33a	MPH00170A	F12	MIMAT0004773	hsa-miR-500	MPH00245A
C01	MIMAT0004702	hsa-miR-339-3p	MPH01266A	G01	MIMAT0002848	hsa-miR-518c	MPH01322A
C02	MIMAT0004958	hsa-miR-301b	MPH00145A	G02	MIMAT0000416	hsa-miR-1	MPH00019A
C03	MIMAT0002819	hsa-miR-193b	MPH01208A	G03	MIMAT0004801	hsa-miR-590-3p	MPH01351A
C04	MIMAT0003322	hsa-miR-652	MPH00408A	G04	MIMAT0000724	hsa-miR-372	MPH00187A
C05	MIMAT0001413	hsa-miR-20b	MPH00105A	G05	MIMAT0000759	hsa-miR-148b	MPH01183A
C06	MIMAT0000430	hsa-miR-138	MPH00039A	G06	MIMAT0000258	hsa-miR-181c	MPH00066A
C07	MIMAT0000232	hsa-miR-199a-3p	MPH01212A	G07	MIMAT0000457	hsa-miR-188-5p	MPH00074A
C08	MIMAT0000261	hsa-miR-183	MPH00069A	G08	MIMAT0000450	hsa-miR-149	MPH00053A
C09	MIMAT0000447	hsa-miR-134	MPH00034A	G09	MIMAT0002171	hsa-miR-410	MPH00205A
C10	MIMAT0000426	hsa-miR-132	MPH01167A	G10	MIMAT0001625	hsa-miR-431	MPH00214A
C11	MIMAT0000243	hsa-miR-148a	MPH01182A	G11	MIMAT0000760	hsa-miR-331-3p	MPH01262A
C12	MIMAT0000275	hsa-miR-218	MPH00115A	G12	MIMAT0004806	hsa-miR-548c-5p	MPH00303A
D01	MIMAT0000710	hsa-miR-365	MPH00182A	H01	MIMAT0002177	hsa-miR-486-5p	MPH00230A
D02	MIMAT0000764	hsa-miR-339-5p	MPH00169A	H02	MIMAT0000686	hsa-miR-34c-5p	MPH00178A
D03	MIMAT0002820	hsa-miR-497	MPH00242A	H03	MIMAT0004810	hsa-miR-629	MPH00385A
D04	MIMAT0004614	hsa-miR-193a-5p	MPH00081A	H04	MIMAT0001635	hsa-miR-452	MPH00223A
D05	MIMAT0003294	hsa-miR-625	MPH00381A	H05	N/A	hsa-SNORD48	MPH01657A
D06	MIMAT0004763	hsa-miR-488	MPH01301A	H06	N/A	hsa-SNORD47	MPH01660A
D07	MIMAT0000242	hsa-miR-129-5p	MPH00028A	H07	N/A	hsa-SNORD44	MPH01658A
D08	MIMAT0000772	hsa-miR-345	MPH00174A	H08	N/A	hsa-U6	MPH01653A
D09	MIMAT0002844	hsa-miR-518b	MPH00266A	H09	N/A	miRTC	N/A
D10	MIMAT0002852	hsa-miR-517a	MPH01318A	H10	N/A	miRTC	N/A
D11	MIMAT0003339	hsa-miR-421	MPH00208A	H11	N/A	PPC	N/A
D12	MIMAT0000227	hsa-miR-197	MPH00087A	H12	N/A	PPC	N/A

Gene Table: MAH-003

Position	mirBASE v10.0 Accession Number	miRNA ID	RT ² miRNA qPCR Assay Cat. #	Position	mirBASE v10.0 Accession Number	miRNA ID	RT ² miRNA qPCR Assay Cat. #
A01	MIMAT0004795	hsa-miR-574-5p	MPH00331A	E01	MIMAT0002880	hsa-miR-508-3p	MPH01311A
A02	MIMAT0000454	hsa-miR-184	MPH00070A	E02	MIMAT0001545	hsa-miR-450a	MPH00220A
A03	MIMAT0003338	hsa-miR-660	MPH00416A	E03	MIMAT0004604	hsa-miR-127-5p	MPH00025A
A04	MIMAT0004695	hsa-miR-337-5p	MPH00167A	E04	MIMAT0002840	hsa-miR-523	MPH01333A
A05	MIMAT0002809	hsa-miR-146b-5p	MPH00048A	E05	MIMAT0004799	hsa-miR-589	MPH00346A
A06	MIMAT0002814	hsa-miR-432	MPH00215A	E06	MIMAT0001621	hsa-miR-369-5p	MPH00184A
A07	MIMAT0000722	hsa-miR-370	MPH00185A	E07	MIMAT0003283	hsa-miR-615-3p	MPH01353A
A08	MIMAT0000094	hsa-miR-95	MPH00478A	E08	MIMAT0002862	hsa-miR-527	MPH00292A
A09	MIMAT0000462	hsa-miR-206	MPH00101A	E09	MIMAT0004796	hsa-miR-576-3p	MPH01348A
A10	MIMAT0004953	hsa-miR-873	MPH00441A	E10	MIMAT0000754	hsa-miR-337-3p	MPH01264A
A11	MIMAT0004780	hsa-miR-532-3p	MPH01337A	E11	MIMAT0004683	hsa-miR-362-3p	MPH01275A
A12	MIMAT0004552	hsa-miR-139-3p	MPH01173A	E12	MIMAT0002827	hsa-miR-515-3p	MPH01315A
B01	MIMAT0002849	hsa-miR-524-5p	MPH00288A	F01	MIMAT0002863	hsa-miR-518a-3p	MPH01321A
B02	MIMAT0004613	hsa-miR-188-3p	MPH01202A	F02	MIMAT0004947	hsa-miR-885-5p	MPH00446A
B03	MIMAT0003258	hsa-miR-590-5p	MPH00347A	F03	MIMAT0001627	hsa-miR-433	MPH00216A
B04	MIMAT0003150	hsa-miR-455-5p	MPH00226A	F04	MIMAT0002823	hsa-miR-512-3p	MPH01313A
B05	MIMAT0004700	hsa-miR-331-5p	MPH00165A	F05	MIMAT0004812	hsa-miR-548d-5p	MPH00304A
B06	MIMAT0001536	hsa-miR-429	MPH00213A	F06	MIMAT0000770	hsa-miR-133b	MPH00033A
B07	MIMAT0000705	hsa-miR-362-5p	MPH00180A	F07	MIMAT0000690	hsa-miR-296-5p	MPH00135A
B08	MIMAT0002839	hsa-miR-525-3p	MPH01335A	F08	MIMAT0000723	hsa-miR-371-3p	MPH01279A
B09	MIMAT0003161	hsa-miR-493	MPH01304A	F09	MIMAT0000751	hsa-miR-330-3p	MPH01261A
B10	MIMAT0002835	hsa-miR-526b	MPH00291A	F10	MIMAT0001629	hsa-miR-329	MPH00163A
B11	MIMAT0000458	hsa-miR-190	MPH00077A	F11	MIMAT0002817	hsa-miR-495	MPH00240A
B12	MIMAT0002826	hsa-miR-515-5p	MPH00261A	F12	MIMAT0004909	hsa-miR-450b-5p	MPH00221A
C01	MIMAT0004952	hsa-miR-665	MPH00420A	G01	MIMAT0003879	hsa-miR-758	MPH00429A
C02	MIMAT0000251	hsa-miR-147	MPH00049A	G02	MIMAT0002868	hsa-miR-522	MPH01332A
C03	MIMAT0004682	hsa-miR-361-3p	MPH01274A	G03	MIMAT0000755	hsa-miR-323-3p	MPH01259A
C04	MIMAT0002883	hsa-miR-514	MPH00260A	G04	MIMAT0004921	hsa-miR-889	MPH00450A
C05	MIMAT0002172	hsa-miR-376b	MPH00194A	G05	MIMAT0004960	hsa-miR-208b	MPH00103A
C06	MIMAT0003247	hsa-miR-582-5p	MPH00339A	G06	MIMAT0002178	hsa-miR-487a	MPH00231A
C07	MIMAT0000262	hsa-miR-187	MPH01201A	G07	MIMAT0003297	hsa-miR-628-3p	MPH01357A
C08	MIMAT0003880	hsa-miR-671-5p	MPH00422A	G08	MIMAT0004809	hsa-miR-628-5p	MPH00384A
C09	MIMAT0003180	hsa-miR-487b	MPH00232A	G09	MIMAT0004819	hsa-miR-671-3p	MPH01360A
C10	MIMAT0003321	hsa-miR-651	MPH00407A	G10	MIMAT0002176	hsa-miR-485-3p	MPH01299A
C11	MIMAT0000276	hsa-miR-219-5p	MPH00116A	G11	MIMAT0003323	hsa-miR-548d-3p	MPH01343A
C12	MIMAT0000687	hsa-miR-299-3p	MPH01243A	G12	MIMAT0000738	hsa-miR-383	MPH00202A
D01	MIMAT0001541	hsa-miR-449a	MPH00218A	H01	MIMAT0003254	hsa-miR-548b-3p	MPH01341A
D02	MIMAT0000721	hsa-miR-369-3p	MPH01278A	H02	MIMAT0000272	hsa-miR-215	MPH00111A
D03	MIMAT0004675	hsa-miR-219-2-3p	MPH01227A	H03	MIMAT0002855	hsa-miR-520d-5p	MPH00280A
D04	MIMAT0000102	hsa-miR-105	MPH00013A	H04	MIMAT0002856	hsa-miR-520d-3p	MPH01331A
D05	MIMAT0003233	hsa-miR-551b	MPH01345A	H05	N/A	hsa-SNORD48	MPH01657A
D06	MIMAT0003340	hsa-miR-542-5p	MPH00296A	H06	N/A	hsa-SNORD47	MPH01660A
D07	MIMAT0003389	hsa-miR-542-3p	MPH01339A	H07	N/A	hsa-SNORD44	MPH01658A
D08	MIMAT0003287	hsa-miR-618	MPH00374A	H08	N/A	hsa-U6	MPH01653A
D09	MIMAT0004775	hsa-miR-502-3p	MPH01309A	H09	N/A	miRTC	N/A
D10	MIMAT0002869	hsa-miR-519a	MPH00272A	H10	N/A	miRTC	N/A
D11	MIMAT0002821	hsa-miR-181d	MPH00067A	H11	N/A	PPC	N/A
D12	MIMAT0000737	hsa-miR-382	MPH00201A	H12	N/A	PPC	N/A

Gene Table: MAH-004

Position	mirBASE v10.0 Accession Number	miRNA ID	RT ² miRNA qPCR Assay Cat. #	Position	mirBASE v10.0 Accession Number	miRNA ID	RT ² miRNA qPCR Assay Cat. #
A01	MIMAT0003249	hsa-miR-584	MPH00341A	E01	MIMAT0002808	hsa-miR-511	MPH00257A
A02	MIMAT0002877	hsa-miR-513-5p	MPH00259A	E02	MIMAT0003263	hsa-miR-595	MPH00351A
A03	MIMAT0003244	hsa-miR-579	MPH00336A	E03	MIMAT0003281	hsa-miR-613	MPH00369A
A04	MIMAT0002824	hsa-miR-498	MPH00243A	E04	MIMAT0004185	hsa-miR-802	MPH00438A
A05	MIMAT0004800	hsa-miR-550	MPH00306A	E05	MIMAT0003248	hsa-miR-583	MPH00340A
A06	MIMAT0003220	hsa-miR-556-5p	MPH00313A	E06	MIMAT0001075	hsa-miR-384	MPH00203A
A07	MIMAT0000735	hsa-miR-380	MPH01286A	E07	MIMAT0003279	hsa-miR-611	MPH00367A
A08	MIMAT0004765	hsa-miR-491-3p	MPH01303A	E08	MIMAT0003240	hsa-miR-575	MPH00332A
A09	MIMAT0003241	hsa-miR-576-5p	MPH00333A	E09	MIMAT0003336	hsa-miR-658	MPH00414A
A10	MIMAT0003265	hsa-miR-597	MPH00353A	E10	MIMAT0003307	hsa-miR-637	MPH00393A
A11	MIMAT0002875	hsa-miR-504	MPH00249A	E11	MIMAT0003245	hsa-miR-580	MPH00337A
A12	MIMAT0002858	hsa-miR-520g	MPH00283A	E12	MIMAT0003228	hsa-miR-564	MPH00321A
B01	MIMAT0003312	hsa-miR-642	MPH00398A	F01	MIMAT0003278	hsa-miR-610	MPH00366A
B02	MIMAT0004922	hsa-miR-875-5p	MPH00443A	F02	MIMAT0003303	hsa-miR-633	MPH00389A
B03	MIMAT0004923	hsa-miR-875-3p	MPH01368A	F03	MIMAT0003226	hsa-miR-562	MPH00319A
B04	MIMAT0002872	hsa-miR-501-5p	MPH00246A	F04	MIMAT0003282	hsa-miR-614	MPH00370A
B05	MIMAT0004774	hsa-miR-501-3p	MPH01308A	F05	MIMAT0000771	hsa-miR-325	MPH00160A
B06	MIMAT0003164	hsa-miR-544	MPH00298A	F06	MIMAT0003268	hsa-miR-600	MPH00356A
B07	MIMAT0004803	hsa-miR-548a-5p	MPH00301A	F07	MIMAT0003299	hsa-miR-630	MPH00386A
B08	MIMAT0003251	hsa-miR-548a-3p	MPH00300A	F08	MIMAT0004901	hsa-miR-298	MPH00137A
B09	MIMAT0001339	hsa-miR-422a	MPH00209A	F09	MIMAT0003291	hsa-miR-622	MPH00378A
B10	MIMAT0002805	hsa-miR-489	MPH00234A	F10	MIMAT0003295	hsa-miR-626	MPH00382A
B11	MIMAT0003306	hsa-miR-636	MPH00392A	F11	MIMAT0003337	hsa-miR-659	MPH00415A
B12	MIMAT0000228	hsa-miR-198	MPH00088A	F12	MIMAT0003250	hsa-miR-585	MPH00342A
C01	MIMAT0002818	hsa-miR-496	MPH00241A	G01	MIMAT0003309	hsa-miR-639	MPH00395A
C02	MIMAT0002864	hsa-miR-518d-3p	MPH01323A	G02	MIMAT0004915	hsa-miR-220c	MPH00120A
C03	MIMAT0002850	hsa-miR-524-3p	MPH01334A	G03	MIMAT0000268	hsa-miR-211	MPH00108A
C04	MIMAT0005454	hsa-miR-519b-5p	MPH00273A	G04	MIMAT0003292	hsa-miR-623	MPH00379A
C05	MIMAT0002832	hsa-miR-519c-3p	MPH01327A	G05	MIMAT0003311	hsa-miR-641	MPH00397A
C06	MIMAT0002860	hsa-miR-516a-3p	MPH01316A	G06	MIMAT0003221	hsa-miR-557	MPH00314A
C07	MIMAT0002873	hsa-miR-502-5p	MPH00247A	G07	MIMAT0000269	hsa-miR-212	MPH00109A
C08	MIMAT0002822	hsa-miR-512-5p	MPH00258A	G08	MIMAT0003302	hsa-miR-632	MPH00388A
C09	MIMAT0003230	hsa-miR-566	MPH00323A	G09	MIMAT0003231	hsa-miR-567	MPH00324A
C10	MIMAT0003285	hsa-miR-548c-3p	MPH01342A	G10	MIMAT0003273	hsa-miR-605	MPH00361A
C11	MIMAT0003219	hsa-miR-555	MPH00312A	G11	MIMAT0000241	hsa-miR-208	MPH00102A
C12	MIMAT0002812	hsa-miR-492	MPH00237A	G12	MIMAT0001532	hsa-miR-448	MPH00217A
D01	MIMAT0003308	hsa-miR-638	MPH00394A	H01	MIMAT0003290	hsa-miR-621	MPH00377A
D02	MIMAT0002170	hsa-miR-412	MPH00207A	H02	MIMAT0003318	hsa-miR-648	MPH00404A
D03	MIMAT0003255	hsa-miR-588	MPH00345A	H03	MIMAT0004929	hsa-miR-190b	MPH00078A
D04	MIMAT0003320	hsa-miR-650	MPH00406A	H04	MIMAT0003253	hsa-miR-587	MPH00344A
D05	MIMAT0003288	hsa-miR-619	MPH00375A	H05	N/A	hsa-SNORD48	MPH01657A
D06	MIMAT0004954	hsa-miR-543	MPH00297A	H06	N/A	hsa-SNORD47	MPH01660A
D07	MIMAT0003313	hsa-miR-643	MPH00399A	H07	N/A	hsa-SNORD44	MPH01658A
D08	MIMAT0003271	hsa-miR-603	MPH00359A	H08	N/A	hsa-U6	MPH01653A
D09	MIMAT0003270	hsa-miR-602	MPH00358A	H09	N/A	miRTC	N/A
D10	MIMAT0003227	hsa-miR-563	MPH00320A	H10	N/A	miRTC	N/A
D11	MIMAT0003222	hsa-miR-558	MPH00315A	H11	N/A	PPC	N/A
D12	MIMAT0003335	hsa-miR-657	MPH00413A	H12	N/A	PPC	N/A